



hacare.com/en/secteur/1/environnement/3/lieu/18/famille/71/collection/298

CHEMODIALYSIS CHAIR-BEDS PICK-UP ARMRESTS - AMB7 UPHOLSTERY - FCD (CLAVIA)


Day care unit - Chemotherapy

Product ref: AMB7-S06P14-S400



HIGHLIGHTS

- Optimize the treatment area
- Double columns and high stability
- Customization of coatings
- Ultra-comfort in a single step
- Two widths to suit every body type
- Facilitate maintenance
- The QR Code for accessing all information




Combating noise pollution

Anti-vibration support, **high-performance shock absorber**: for **smooth, quiet, and comfortable transport**, even when crossing thresholds or exiting elevators. **Noise pollution is significantly reduced.**



Combating musculoskeletal disorders

Ergonomic handles, double casters, and a **fifth synchronized swivel wheel** for **optimal maneuverability**. Designed to reduce musculoskeletal disorders and adapt to all body types.



Make cleaning easier

Disassembly and cover removal (via buttonholes) **without tools** for quick, **thorough cleaning, perfect hygiene**, and real comfort for nursing staff.



Ergonomics and comfort during sampling or treatment

Designed to fit all body types, **Body Flex armrests** are **multi-directional, adjustable**, and **can be moved 7.5 cm** for optimal arm positioning.

Performance

• Length sleeping surface	202 cm	CE marking	CE
• Width sleeping surface	65 cm		
• Variable height minimum	Sleeping surface without mattress: 47 cm		
• Variable height maximum	Sleeping surface without mattress: 85 cm		
• Foam HR density	35 kg/m ³		
• Backrest	0 to 78°		
• Legrest	0 to 54°		
• Trendelenburg	15°		
• Maximum patient weight	185 kg		









Improve according to your needs

 17 Options

 10 Accessories

 20 Color options

Complementary products



Advisors close to you

Export HAcare
export@hacare.com
+33 (0)4 74 08 71 71