



hacare.com/en/secteur/1/environnement/3/lieu/59/famille/71/collection/326


CHEMODIALYSIS CHAIR-BEDS SAMPLING ARMRESTS - AMB7 MATTRESS - FCD (ERGOLYS)

Day care unit - Dialysis
Product ref: AMBE-S07P14-S400



HIGHLIGHTS

- Long-lasting comfort and optimal ergonomics
- Double columns and high stability
- Customized coatings
- Optimizing the care surface
- Easy cleaning
- Easy maintenance
- The QR Code to access any information



Fighting noise pollution

Anti-vibration support, **high-performance damper**: for smooth, **quiet** and **comfortable** stretching, even when crossing thresholds or exiting elevators. **Noise pollution is significantly reduced**



Ultra-comfortable armrests for picking or care

Designed to adjust to all body types, the **Body Flex armrests** are **multi-directional**, **orientable** and **translatable** by 7.5 cm for optimal arm positioning.



Fighting MSDs

Ergonomic handles, double-roller wheels and **5th directional and synchronized wheel** for optimal **maneuverability**. Designed to limit RSI and adapt to all morphologies.



Ultra-comfort in a single gesture

Thanks to the "**Zero gravity position**", **recognized** for its ability to **distribute body pressure harmoniously**, the patient's weight is distributed evenly over the entire mattress. This posture **effectively reduces pressure points** and provides a pleasant sensation of lightness and well-being.
Optimized comfort, relieved tension,

Performance

• Length sleeping surface	213 cm	Maximum patient weight	200 kg
• Width sleeping surface	73 cm	CE marking	CE
• Variable height minimum	Sleeping surface without mattress: 47 cm		
• Variable height maximum	Sleeping surface without mattress: 85 cm		
• Foam HR density	37 kg/m ³		
• Viscoelastic foam density	40 kg/m ³		
• Backrest	0 to 78°		
• Legrest	0 to 54°		
• Trendelenburg	15°		



Improve according to your needs

 **16 Options**

 **8 Accessories**

 **20 Color options**

Complementary products



Advisors close to you

Export HAcare
export@hacare.com
+33 (0)4 74 08 71 71